

SURGICAL PRIMADO SERIES CUTTING ACCESSORIES



Create it 

Primado Series

Cutting Accessories



Attachment/Bur Compatibility Table

For Slim Attachment 300

Slim Tube		Bur		Ring Mark
Order Code	Type	Order Code*		
P300-1T110 P300-1T110-C P300-1T110-PB	110	PDS-1□□110-□□		
P300-1T130 P300-1T130-C P300-1T130-PB	130	PDS-1□□130-□□		
P300-1T170	170	PDS-1□□170-□□		
P300-1T240	240	PDS-1□□240-□□		

For Standard Attachment 300

Attachment		Bur		Ring Mark
Order Code	Type	Order Code*		
P300-2S10 P300-2S20 P300-2S40 P300-2S60 P300-2S80	Universal Length	PDS-2□□UL-□□		
P300-2A10 P300-2A20 P300-2A40 P300-2A60 P300-2A80				

For Super Slim Handpiece 200/Super Slim Handpiece

Attachment		Bur		Ring Mark
Order Code	Type	Order Code*		
P200-1A135	L	PDS-1□□135-□□		
P200-1A155 P200-1A155-C	EL	PDS-1□□155-□□		
P200-1A175 P200-1A175-C	SL	PDS-1□□175-□□		
P200-RA330 P200-RA330-L	SPL	PDS-1□□330-□□	No Ring Mark	

For Standard Handpiece 200

Attachment		Bur		Ring Mark
Order Code	Type	Order Code*		
P200-2SMS P200-2AMS	MS	PDS-2□□MS-□□		
P200-2SES P200-2AES	ES	PDS-2□□ES-□□		
P200-2SS P200-2AS	S	PDS-2□□S-□□		
P200-2SM P200-2AM	M	PDS-2□□M-□□		
P200-2SL P200-2AL	L	PDS-2□□L-□□		

For High Speed Handpiece 200

Attachment		Bur		Ring Mark
Order Code	Type	Order Code*		
P200-3GAS	S	PDS-3□□S-□□		
P200-3GAM	M	PDS-3□□M-□□		
P200-3GAL	L	PDS-3□□L-□□		

*□□ represents an alphanumeric character of each Order Code.

Contents

The color index indicates the compatible cutting accessories with the attachments.

Slim Attachment 300	P300-1AHS P300-1AHA	P300-1T110 P300-1T110-C P300-1T110-PB	P300-1T130 P300-1T130-C P300-1T130-PB	P300-1T170 P300-1T240
Standard Attachment 300	P300-2S10 P300-2S20 P300-2S40	P300-2S60 P300-2S80	P300-2A10 P300-2A20 P300-2A40	P300-2A60 P300-2A80
Slim Handpiece 200	P200-1A135 P200-1A155 P200-1A175	P200-1A155-C P200-1A175-C		
Super Slim Handpiece 200 Super Slim Handpiece	P200-RA330 P200-RA330-L	PD-RA330 PD-RA330-L		
Standard Handpiece 200 Standard Handpiece	P200-2SMS P200-2SES P200-2SS	P200-2SM P200-2SL	P200-2AMS P200-2AES P200-2AS	P200-2AM P200-2AL PD-2SES PD-2SS PD-2SM PD-2AS PD-2AM PD-2AL
Drill Handpiece 200 Drill Handpiece	P200-2SD	PD-2SD PD-2AD		
High Speed Handpiece 200	P200-3GAS P200-3GAM P200-3GAL			
Craniotome Handpiece 200 Craniotome Handpiece	P200-CRA P200-CRD-P P200-CRD-M P200-CRD-L	P200-RCRD-P P200-RCRD-M P200-RCRD-L	PD-CRA PD-CRD-P PD-CRD-M PD-CRD-L	
Metal Cutting Handpiece 200	P200-3MC			
Reciprocating Saw Attachment	P200-REC	PD-SA-R		
Sagittal Saw Attachment	P200-SAG	PD-SA-S		
Oscillating Saw Attachment	P200-OSC			

Rotary Instruments For Slim Attachment 300

R Diamond Bur



Head Size (mm)	Slim Tube			
	110	130	170	240
2.0	PDS-1RD110-20	PDS-1RD130-20	—	—
3.0	PDS-1RD110-30	PDS-1RD130-30	—	PDS-1RD240-30
3.5	—	—	PDS-1RD170-35	—
4.0	PDS-1RD110-40	PDS-1RD130-40	PDS-1RD170-40	PDS-1RD240-40
4.5	—	—	—	PDS-1RD240-45
5.0	—	—	PDS-1RD170-50	PDS-1RD240-50

R Coarse Diamond Bur



Head Size (mm)	Slim Tube			
	110	130	170	240
2.0	PDS-1RCD110-20	PDS-1RCD130-20	—	—
3.0	PDS-1RCD110-30	PDS-1RCD130-30	—	PDS-1RCD240-30
3.5	—	—	PDS-1RCD170-35	—
4.0	PDS-1RCD110-40	PDS-1RCD130-40	PDS-1RCD170-40	PDS-1RCD240-40
4.5	—	—	—	PDS-1RCD240-45
5.0	PDS-1RCD110-50	PDS-1RCD130-50	PDS-1RCD170-50	PDS-1RCD240-50

MH Diamond Bur



Head Size (mm)	Slim Tube	
	110	130
2.0	PDS-1MHD110-20	PDS-1MHD130-20
3.0	PDS-1MHD110-30	PDS-1MHD130-30
4.0	PDS-1MHD110-40	PDS-1MHD130-40

RT Diamond Bur



Head Size (mm)	Slim Tube	
	110	130
4.5	PDS-1RTD110-45	PDS-1RTD130-45

MH Coarse Diamond Bur



Head Size (mm)	Slim Tube	
	110	130
3.0	PDS-1MHCD110-30	PDS-1MHCD130-30
4.0	PDS-1MHCD110-40	PDS-1MHCD130-40

RT Coarse Diamond Bur



Head Size (mm)	Slim Tube	
	110	130
4.5	PDS-1RTCD110-45	PDS-1RTCD130-45

Rotary Instruments For Slim Attachment 300

R Fluted Bur



Head Size (mm)	Slim Tube			
	110	130	170	240
3.0	PDS-1RF110-30	PDS-1RF130-30	—	—
3.5	—	—	PDS-1RF170-35	—
4.0	PDS-1RF110-40	PDS-1RF130-40	PDS-1RF170-40	PDS-1RF240-40

MH Fluted Bur



Head Size (mm)	Slim Tube	
	110	130
2.5	PDS-1MHF110-25	PDS-1MHF130-25
3.0	PDS-1MHF110-30	PDS-1MHF130-30

Rotary Instruments For Standard Attachment 300

R Diamond Bur UL



Head Size (mm)	Attachment
	UL
0.6	PDS-2RDUL-06
2.0	PDS-2RDUL-20
2.5	PDS-2RDUL-25
3.0	PDS-2RDUL-30
3.5	PDS-2RDUL-35
4.0	PDS-2RDUL-40
5.0	PDS-2RDUL-50
6.0	PDS-2RDUL-60

R Coarse Diamond Bur UL



Head Size (mm)	Attachment
	UL
1.0	PDS-2RCDUL-10
1.5	PDS-2RCDUL-15
2.0	PDS-2RCDUL-20
2.5	PDS-2RCDUL-25
3.0	PDS-2RCDUL-30
4.0	PDS-2RCDUL-40
5.0	PDS-2RCDUL-50
6.0	PDS-2RCDUL-60
7.0	PDS-2RCDUL-70

R Extra Coarse Diamond Bur UL



Head Size (mm)	Attachment
	UL
3.0	PDS-2RECDUL-30
4.0	PDS-2RECDUL-40
5.0	PDS-2RECDUL-50
6.0	PDS-2RECDUL-60
7.0	PDS-2RECDUL-70

R Super Coarse Diamond Bur UL



Head Size (mm)	Attachment
	UL
3.0	PDS-2RSCDUL-30
4.0	PDS-2RSCDUL-40
5.0	PDS-2RSCDUL-50
6.0	PDS-2RSCDUL-60

MH Coarse Diamond Bur UL



Head Size (mm)	Attachment
	UL
2.0	PDS-2MHCDUL-20
3.0	PDS-2MHCDUL-30
4.0	PDS-2MHCDUL-40

Barrel Diamond Bur UL



Head Size (mm)	Attachment
	UL
3.0	PDS-2BCDUL-30
5.0	PDS-2BCDUL-50

Rotary Instruments For Standard Attachment 300

R Fluted Bur UL



Head Size (mm)	Attachment
	UL
1.0	PDS-2RFUL-10
1.5	PDS-2RFUL-15
2.0	PDS-2RFUL-20
2.5	PDS-2RFUL-25
3.0	PDS-2RFUL-30
3.5	PDS-2RFUL-35
4.0	PDS-2RFUL-40
5.0	PDS-2RFUL-50
6.0	PDS-2RFUL-60
7.0	PDS-2RFUL-70
8.8	PDS-2RFUL-88

QuadEdge Bur UL*



Head Size (mm)	Attachment
	UL
3.0	PDS-2QEUL-30
4.0	PDS-2QEUL-40
5.0	PDS-2QEUL-50
6.0	PDS-2QEUL-60

MH Fluted Bur UL



Head Size (mm)	Attachment
	UL
1.8	PDS-2MHFUL-18
2.2	PDS-2MHFUL-22
3.0	PDS-2MHFUL-30

Egg Fluted Bur UL



Head Size (mm)	Attachment
	UL
4.0	PDS-2EFUL-40
6.0	PDS-2EFUL-60

Cylinder Bur UL



Head Size (mm)	Attachment
	UL
5.0	PDS-2CYUL-50
6.0	PDS-2CYUL-60

Acorn Bur UL*



Head Size (mm)	Attachment
	UL
5.0	PDS-2AUL-50
6.0	PDS-2AUL-60
7.5	PDS-2AUL-75

Router Bur UL*



Head Size (mm)	Attachment
	UL
1.3	PDS-2RUL-13

Wirepass Drill UL*



Head Size (mm)	Attachment
	UL
1.0	PDS-2WUL-10
1.2	PDS-2WUL-12
1.5	PDS-2WUL-15

*The Exposure Length of these burs are not adjustable.

Rotary Instruments For Slim Handpiece 200

Diamond Bur



Head Size (mm)	Attachment		
	L	EL	SL
2.0	PDS-1D135-20	PDS-1D155-20	PDS-1D175-20
3.0	PDS-1D135-30	PDS-1D155-30	PDS-1D175-30
4.0	PDS-1D135-40	PDS-1D155-40	PDS-1D175-40

Coarse Diamond Bur



Head Size (mm)	Attachment				
	L	EL	EL Extended	SL	SL Extended
2.0	PDS-1CD135-20	PDS-1CD155-20	PDS-1CD155-20E	PDS-1CD175-20	PDS-1CD175-20E
3.0	PDS-1CD135-30	PDS-1CD155-30	PDS-1CD155-30E	PDS-1CD175-30	PDS-1CD175-30E
4.0	PDS-1CD135-40	PDS-1CD155-40	—	PDS-1CD175-40	—

Extended The burs can be extended by 6 mm from the original length.

Neuro Diamond Bur



Head Size (mm)	Attachment		
	L	EL	SL
2.0	PDS-1ND135-20	PDS-1ND155-20	PDS-1ND175-20
3.0	PDS-1ND135-30	PDS-1ND155-30	PDS-1ND175-30
4.0	PDS-1ND135-40	PDS-1ND155-40	PDS-1ND175-40

Neuro Coarse Diamond Bur



Head Size (mm)	Attachment		
	L	EL	SL
2.0	PDS-1NCD135-20	PDS-1NCD155-20	PDS-1NCD175-20
3.0	PDS-1NCD135-30	PDS-1NCD155-30	PDS-1NCD175-30
4.0	PDS-1NCD135-40	PDS-1NCD155-40	PDS-1NCD175-40

Round Fluted Bur



Head Size (mm)	Attachment	
	EL	SL
2.0	PDS-1C155-20	PDS-1C175-20
3.0	PDS-1C155-30	PDS-1C175-30
4.0	PDS-1C155-40	PDS-1C175-40

MH Cutting Bur



Head Size (mm)	Attachment	
	EL	SL
2.5	PDS-1MHC155-25	PDS-1MHC175-25
3.0	PDS-1MHC155-30	PDS-1MHC175-30

For Super Slim Handpiece 200 Super Slim Handpiece

Diamond Bur



Head Size (mm)	Attachment
	SPL
2.0	PDS-1D330-20
2.5	PDS-1D330-25
3.0	PDS-1D330-30
3.5	PDS-1D330-35

Coarse Diamond Bur



Head Size (mm)	Attachment
	SPL
2.0	PDS-1CD330-20
2.5	PDS-1CD330-25
3.0	PDS-1CD330-30
3.5	PDS-1CD330-35

Rotary Instruments For Standard Handpiece 200/Standard Handpiece

Diamond Bur



Head Size (mm)	Attachment				
	MS	ES	S	M	L
0.6	PDS-2DMS-06	PDS-2DES-06	PDS-2DS-06	—	—
1.0	PDS-2DMS-10	PDS-2DES-10	PDS-2DS-10	PDS-2DM-10	—
1.5	PDS-2DMS-15	PDS-2DES-15	PDS-2DS-15	—	—
2.0	PDS-2DMS-20	PDS-2DES-20	PDS-2DS-20	PDS-2DM-20	PDS-2DL-20
3.0	PDS-2DMS-30	PDS-2DES-30	PDS-2DS-30	PDS-2DM-30	PDS-2DL-30
4.0	PDS-2DMS-40	PDS-2DES-40	PDS-2DS-40	PDS-2DM-40	PDS-2DL-40
5.0	PDS-2DMS-50	PDS-2DES-50	PDS-2DS-50	PDS-2DM-50	PDS-2DL-50
6.0	PDS-2DMS-60	PDS-2DES-60	PDS-2DS-60	PDS-2DM-60	PDS-2DL-60
7.0	PDS-2DMS-70	PDS-2DES-70	PDS-2DS-70	PDS-2DM-70	PDS-2DL-70

Coarse Diamond Bur



Head Size (mm)	Attachment				
	MS	ES	S	M	L
1.0	PDS-2CDMS-10	PDS-2CDES-10	PDS-2CDS-10	—	—
1.5	PDS-2CDMS-15	PDS-2CDES-15	PDS-2CDS-15	—	—
2.0	PDS-2CDMS-20	PDS-2CDES-20	PDS-2CDS-20	PDS-2CDM-20	PDS-2CDL-20
3.0	PDS-2CDMS-30	PDS-2CDES-30	PDS-2CDS-30	PDS-2CDM-30	PDS-2CDL-30
4.0	PDS-2CDMS-40	PDS-2CDES-40	PDS-2CDS-40	PDS-2CDM-40	PDS-2CDL-40
5.0	PDS-2CDMS-50	PDS-2CDES-50	PDS-2CDS-50	PDS-2CDM-50	PDS-2CDL-50
6.0	PDS-2CDMS-60	PDS-2CDES-60	PDS-2CDS-60	PDS-2CDM-60	PDS-2CDL-60
7.0	PDS-2CDMS-70	PDS-2CDES-70	PDS-2CDS-70	PDS-2CDM-70	PDS-2CDL-70

Coarse Diamond Bur/Extended



Head Size (mm)	Attachment			
	ES Extended	S Extended	M Extended	L Extended
1.0	PDS-2CDES-10E	—	—	—
1.5	PDS-2CDES-15E	—	—	—
2.0	—	PDS-2CDS-20E	PDS-2CDM-20E	PDS-2CDL-20E
3.0	—	PDS-2CDS-30E	PDS-2CDM-30E	PDS-2CDL-30E
4.0	—	PDS-2CDS-40E	PDS-2CDM-40E	PDS-2CDL-40E
5.0	—	PDS-2CDS-50E	PDS-2CDM-50E	PDS-2CDL-50E

Extended The burs can be extended by 10 mm from the original length.

Extra Coarse Diamond Bur



Head Size (mm)	Attachment				
	MS	ES	S	M	L
3.0	PDS-2ECDMS-30	PDS-2ECDES-30	PDS-2ECDS-30	PDS-2ECDM-30	PDS-2ECDL-30
4.0	PDS-2ECDMS-40	PDS-2ECDES-40	PDS-2ECDS-40	PDS-2ECDM-40	PDS-2ECDL-40
5.0	PDS-2ECDMS-50	PDS-2ECDES-50	PDS-2ECDS-50	PDS-2ECDM-50	PDS-2ECDL-50
6.0	PDS-2ECDMS-60	PDS-2ECDES-60	PDS-2ECDS-60	PDS-2ECDM-60	PDS-2ECDL-60
7.0	PDS-2ECDMS-70	—	—	—	—

Rotary Instruments

For Standard Handpiece 200/Standard Handpiece

Super Coarse Diamond Bur



Head Size (mm)	Attachment				
	MS	ES	S	M	L
3.0	PDS-2SCDMS-30	PDS-2SCDES-30	PDS-2SCDS-30	PDS-2SCDM-30	PDS-2SCDL-30
4.0	PDS-2SCDMS-40	PDS-2SCDES-40	PDS-2SCDS-40	PDS-2SCDM-40	PDS-2SCDL-40
5.0	PDS-2SCDMS-50	PDS-2SCDES-50	PDS-2SCDS-50	PDS-2SCDM-50	PDS-2SCDL-50
6.0	PDS-2SCDMS-60	PDS-2SCDES-60	PDS-2SCDS-60	PDS-2SCDM-60	PDS-2SCDL-60

Barrel Diamond Bur



Head Size (mm)	Attachment	
	L	L <small>Extended</small>
3.0	PDS-2BCDL-30	PDS-2BCDL-30E
5.0	PDS-2BCDL-50	PDS-2BCDL-50E

Extended The burs can be extended by 10 mm from the original length.

Round Fluted Bur



Head Size (mm)	Attachment				
	MS	ES	S	M	L
1.0	PDS-2CMS-10	PDS-2CES-10	PDS-2CS-10	PDS-2CM-10	—
2.0	PDS-2CMS-20	PDS-2CES-20	PDS-2CS-20	PDS-2CM-20	PDS-2CL-20
3.0	PDS-2CMS-30	PDS-2CES-30	PDS-2CS-30	PDS-2CM-30	PDS-2CL-30
4.0	PDS-2CMS-40	PDS-2CES-40	PDS-2CS-40	PDS-2CM-40	PDS-2CL-40
5.0	PDS-2CMS-50	PDS-2CES-50	PDS-2CS-50	PDS-2CM-50	PDS-2CL-50
6.0	PDS-2CMS-60	PDS-2CES-60	PDS-2CS-60	PDS-2CM-60	PDS-2CL-60
7.0	PDS-2CMS-70	PDS-2CES-70	PDS-2CS-70	PDS-2CM-70	PDS-2CL-70

Carbide Bur



Head Size (mm)	Attachment				
	MS	ES	S	M	L
2.0	PDS-2CBMS-20	PDS-2CBES-20	PDS-2CBS-20	PDS-2CBM-20	PDS-2CBL-20
3.0	PDS-2CBMS-30	PDS-2CBES-30	PDS-2CBS-30	PDS-2CBM-30	PDS-2CBL-30
4.0	PDS-2CBMS-40	PDS-2CBES-40	PDS-2CBS-40	PDS-2CBM-40	PDS-2CBL-40
5.0	PDS-2CBMS-50	PDS-2CBES-50	PDS-2CBS-50	PDS-2CBM-50	PDS-2CBL-50
6.0	PDS-2CBMS-60	PDS-2CBES-60	PDS-2CBS-60	PDS-2CBM-60	PDS-2CBL-60
7.0	PDS-2CBMS-70	PDS-2CBES-70	PDS-2CBS-70	PDS-2CBM-70	PDS-2CBL-70

Neuro Bur



Head Size (mm)	Attachment				
	MS	ES	S	M	L
1.8	PDS-2NMS-18	PDS-2NES-18	PDS-2NS-18	PDS-2NM-18	PDS-2NL-18
2.2	PDS-2NMS-22	PDS-2NES-22	PDS-2NS-22	PDS-2NM-22	PDS-2NL-22
3.0	PDS-2NMS-30	PDS-2NES-30	PDS-2NS-30	PDS-2NM-30	PDS-2NL-30

Rotary Instruments

For Standard Handpiece 200/Standard Handpiece

Neuro Carbide Bur



Head Size (mm)	Attachment				
	MS	ES	S	M	L
1.8	PDS-2NCBMS-18	PDS-2NCBES-18	PDS-2NCBS-18	PDS-2NCBM-18	PDS-2NCBL-18
2.2	PDS-2NCBMS-22	PDS-2NCBES-22	PDS-2NCBS-22	PDS-2NCBM-22	PDS-2NCBL-22
3.0	PDS-2NCBMS-30	PDS-2NCBES-30	PDS-2NCBS-30	PDS-2NCBM-30	PDS-2NCBL-30

Egg Bur



Head Size (mm)	Attachment			
	MS	ES	S	L
4.0	PDS-2EMS-40	PDS-2EES-40	PDS-2ES-40	PDS-2EL-40
6.0	PDS-2EMS-60	PDS-2EES-60	PDS-2ES-60	PDS-2EL-60

Barrel Bur



Head Size (mm)	Attachment		
	MS	ES	S
4.0	PDS-2BMS-40	PDS-2BES-40	PDS-2BS-40
6.0	PDS-2BMS-60	PDS-2BES-60	PDS-2BS-60

Pear Bur



Head Size (mm)	Attachment		
	MS	ES	S
6.0	PDS-2PMS-60	PDS-2PES-60	PDS-2PS-60

Cylinder Bur



Head Size (mm)	Attachment
	L
5.0	PDS-2CYL-50
6.0	PDS-2CYL-60

Acorn Bur



Head Size (mm)	Attachment		
	MS	ES	S
5.0	PDS-2AMS-50	PDS-2AES-50	PDS-2AS-50
6.0	PDS-2AMS-60	PDS-2AES-60	PDS-2AS-60
7.5	PDS-2AMS-75	PDS-2AES-75	PDS-2AS-75

Router Bur



Head Size (mm)	Attachment	
	MS	ES
1.3	PDS-2RMS-13	PDS-2RES-13

Lindemann Bur



Blade Length x Head Size (mm)	Attachment		
	MS	ES	S
11x2.1	PDS-2LMS-11	PDS-2LES-11	PDS-2LS-11
26x2.3	PDS-2LMS-26	PDS-2LES-26	PDS-2LS-26

Fissure Bur



Head Size (mm)	Attachment		
	MS	ES	S
1.2	PDS-2FMS-12	PDS-2FES-12	PDS-2FS-12
1.6	PDS-2FMS-16	PDS-2FES-16	PDS-2FS-16
2.1	PDS-2FMS-21	PDS-2FES-21	PDS-2FS-21

Wirepass Drill



Head Size (mm)	Attachment		
	MS	ES	S
1.0	PDS-2WMS-10	PDS-2WES-10	PDS-2WS-10
1.2	PDS-2WMS-12	PDS-2WES-12	PDS-2WS-12
1.5	PDS-2WMS-15	PDS-2WES-15	PDS-2WS-15

Rotary Instruments

For Drill Handpiece 200/Drill Handpiece

Twist Drill Flat Tail



Head Size (mm)	Attachment
	2SD/2AD
1.0	PDS-2TD-10
1.5	PDS-2TD-15
2.0	PDS-2TD-20

For High Speed Handpiece 200

Diamond Bur



Head Size (mm)	Attachment				
	S	M	M Extended	L	L Extended
1.0	PDS-3DS-10	PDS-3DM-10	—	—	—
2.0	PDS-3DS-20	PDS-3DM-20	PDS-3DM-20E	PDS-3DL-20	PDS-3DL-20E
3.0	PDS-3DS-30	PDS-3DM-30	—	PDS-3DL-30	—
4.0	PDS-3DS-40	PDS-3DM-40	—	PDS-3DL-40	—
5.0	PDS-3DS-50	PDS-3DM-50	—	PDS-3DL-50	—
6.0	PDS-3DS-60	PDS-3DM-60	—	PDS-3DL-60	—
7.0	PDS-3DS-70	PDS-3DM-70	—	PDS-3DL-70	—

Extended The burs can be extended by 10 mm from the original length.

Coarse Diamond Bur



Head Size (mm)	Attachment				
	S	M	M Extended	L	L Extended
2.0	PDS-3CDS-20	PDS-3CDM-20	—	PDS-3CDL-20	—
3.0	PDS-3CDS-30	PDS-3CDM-30	PDS-3CDM-30E	PDS-3CDL-30	PDS-3CDL-30E
4.0	PDS-3CDS-40	PDS-3CDM-40	PDS-3CDM-40E	PDS-3CDL-40	PDS-3CDL-40E
5.0	PDS-3CDS-50	PDS-3CDM-50	PDS-3CDM-50E	PDS-3CDL-50	PDS-3CDL-50E
6.0	PDS-3CDS-60	PDS-3CDM-60	—	PDS-3CDL-60	—
7.0	PDS-3CDS-70	PDS-3CDM-70	—	PDS-3CDL-70	—

Extended The burs can be extended by 10 mm from the original length.

Round Fluted Bur



Head Size (mm)	Attachment		
	S	M	L
1.0	PDS-3CS-10	PDS-3CM-10	—
2.0	PDS-3CS-20	PDS-3CM-20	PDS-3CL-20
3.0	PDS-3CS-30	PDS-3CM-30	PDS-3CL-30
4.0	PDS-3CS-40	PDS-3CM-40	PDS-3CL-40
5.0	PDS-3CS-50	PDS-3CM-50	PDS-3CL-50
6.0	PDS-3CS-60	PDS-3CM-60	PDS-3CL-60
7.0	PDS-3CS-70	PDS-3CM-70	PDS-3CL-70

Neuro Bur



Head Size (mm)	Attachment		
	S	M	L
1.8	PDS-3NS-18	PDS-3NM-18	PDS-3NL-18
2.2	PDS-3NS-22	PDS-3NM-22	PDS-3NL-22
3.0	PDS-3NS-30	PDS-3NM-30	PDS-3NL-30

Rotary Instruments

For Craniotome Handpiece 200/ Craniotome Handpiece

Craniotome Handpiece Bur



Attachment		
P	M	L
PDS-CRA-P	PDS-CRA-M	PDS-CRA-L

Router Bur



Attachment		
P	M	L
PDS-ROU-P	PDS-ROU-M	PDS-ROU-L

Waved Spiral Bur



Attachment	
M	L
PDS-WSP-M	PDS-WSP-L

For Metal Cutting Handpiece 200

Metal Cutting Carbide Bur



Head Size (mm)	Attachment
	3MC
1.6	PDS-3MCS-CB16



Head Size (mm)	Attachment
	3MC
3.0	PDS-3MCS-CB30

Metal Cutting Diamond Disc




Head Size (mm)	Attachment
	3MC
25.4	PDS-3MCS-DD254

Saw Blades

Reciprocating Saw Blade

Scale: 1/1



PDS-RE-4023

Blade Length	23.4 mm
Thickness	0.38 mm

PDS-RE-4023E

Blade Length	23.4 mm
Thickness	0.38 mm

Extended The shank can be extended by 10 mm from the original length.

PDS-RE-4028

Blade Length	28.0 mm
Thickness	0.38 mm

PDS-RE-4028E

Blade Length	28.0 mm
Thickness	0.38 mm

Extended The shank can be extended by 10 mm from the original length.

PDS-RE-4114

Blade Length	14.5 mm
Thickness	0.4 mm

PDS-RE-4127

Blade Length	27.0 mm
Thickness	0.4 mm

PDS-RE-4217

Blade Length	17.5 mm
Thickness	0.4 mm

PDS-RE-4223

Blade Length	23.0 mm
Thickness	0.4 mm


PDS-RE-6127

Blade Length	27.0 mm
Thickness	0.6 mm

Rasps

Reciprocating Rasp

Scale: 1/1



PDS-RE-R011

· Straight-Cut

Blade Length	11.0 mm
Thickness	5.5 mm

PDS-RE-R012

· Straight-Cut

Blade Length	13.0 mm
Thickness	7.0 mm

PDS-RE-R111

· Cross-Cut

Blade Length	11.0 mm
Thickness	5.3 mm

PDS-RE-R112


· Cross-Cut

Blade Length	13.0 mm
Thickness	7.0 mm

Saw Blades


Sagittal and Oscillating Saw Blade

Scale: 1/1




PDS-SO-411

Cutting Length	11.5 mm
Width	5.5 mm
Thickness	0.4 mm




PDS-SO-442

Cutting Length	25.0 mm
Width	9.3 mm
Thickness	0.4 mm




PDS-SO-412

Cutting Length	11.5 mm
Width	9.0 mm
Thickness	0.4 mm




PDS-SO-452

Cutting Length	31.0 mm
Width	9.0 mm
Thickness	0.4 mm




PDS-SO-422

Cutting Length	15.0 mm
Width	7.0 mm
Thickness	0.4 mm




PDS-SO-453

Cutting Length	39.0 mm
Width	13.2 mm
Thickness	0.4 mm




PDS-SO-431

Cutting Length	18.0 mm
Width	5.5 mm
Thickness	0.4 mm




PDS-SO-632

Cutting Length	18.0 mm
Width	9.0 mm
Thickness	0.6 mm




PDS-SO-432

Cutting Length	18.0 mm
Width	9.0 mm
Thickness	0.4 mm




PDS-SO-642

Cutting Length	25.0 mm
Width	9.3 mm
Thickness	0.6 mm




PDS-SO-433

Cutting Length	18.5 mm
Width	16.5 mm
Thickness	0.4 mm



PDS-SO-653

Cutting Length	39.0 mm
Width	13.2 mm
Thickness	0.6 mm





PDS-SO-441

Cutting Length	25.0 mm
Width	5.5 mm
Thickness	0.4 mm

Saw Blades



Oscillating Saw Blade

Scale: 1/1



PDS-OSC-422*

Cutting Length	7.0 mm
Width	11.5 mm
Thickness	0.4 mm



PDS-OSC-432*

Cutting Length	12.0 mm
Width	11.5 mm
Thickness	0.4 mm

*PD-OSC-A is required for Oscillating Saw Blades installation.

Index

Order Code	Page
PDS-1C155-20	p.07
PDS-1C155-30	p.07
PDS-1C155-40	p.07
PDS-1C175-20	p.07
PDS-1C175-30	p.07
PDS-1C175-40	p.07
PDS-1CD135-20	p.07
PDS-1CD135-30	p.07
PDS-1CD135-40	p.07
PDS-1CD155-20	p.07
PDS-1CD155-20E	p.07
PDS-1CD155-30	p.07
PDS-1CD155-30E	p.07
PDS-1CD155-40	p.07
PDS-1CD175-20	p.07
PDS-1CD175-20E	p.07
PDS-1CD175-30	p.07
PDS-1CD175-30E	p.07
PDS-1CD175-40	p.07
PDS-1CD330-20	p.07
PDS-1CD330-25	p.07
PDS-1CD330-30	p.07
PDS-1CD330-35	p.07
PDS-1D135-20	p.07
PDS-1D135-30	p.07
PDS-1D135-40	p.07
PDS-1D155-20	p.07
PDS-1D155-30	p.07
PDS-1D155-40	p.07
PDS-1D175-20	p.07
PDS-1D175-30	p.07
PDS-1D175-40	p.07
PDS-1D330-20	p.07
PDS-1D330-25	p.07
PDS-1D330-30	p.07
PDS-1D330-35	p.07
PDS-1MHC155-25	p.07
PDS-1MHC155-30	p.07
PDS-1MHC175-25	p.07
PDS-1MHC175-30	p.07
PDS-1MHCD110-30	p.03
PDS-1MHCD110-40	p.03
PDS-1MHCD130-30	p.03
PDS-1MHCD130-40	p.03
PDS-1MHD110-20	p.03
PDS-1MHD110-30	p.03
PDS-1MHD110-40	p.03
PDS-1MHD130-20	p.03
PDS-1MHD130-30	p.03
PDS-1MHD130-40	p.03

Order Code	Page
PDS-1MHF110-25	p.04
PDS-1MHF110-30	p.04
PDS-1MHF130-25	p.04
PDS-1MHF130-30	p.04
PDS-1NCD135-20	p.07
PDS-1NCD135-30	p.07
PDS-1NCD135-40	p.07
PDS-1NCD155-20	p.07
PDS-1NCD155-30	p.07
PDS-1NCD155-40	p.07
PDS-1NCD175-20	p.07
PDS-1NCD175-30	p.07
PDS-1NCD175-40	p.07
PDS-1ND135-20	p.07
PDS-1ND135-30	p.07
PDS-1ND135-40	p.07
PDS-1ND155-20	p.07
PDS-1ND155-30	p.07
PDS-1ND155-40	p.07
PDS-1ND175-20	p.07
PDS-1ND175-30	p.07
PDS-1ND175-40	p.07
PDS-1RCD110-20	p.03
PDS-1RCD110-30	p.03
PDS-1RCD110-40	p.03
PDS-1RCD110-50	p.03
PDS-1RCD130-20	p.03
PDS-1RCD130-30	p.03
PDS-1RCD130-40	p.03
PDS-1RCD130-50	p.03
PDS-1RCD170-35	p.03
PDS-1RCD170-40	p.03
PDS-1RCD170-50	p.03
PDS-1RCD240-30	p.03
PDS-1RCD240-40	p.03
PDS-1RCD240-45	p.03
PDS-1RCD240-50	p.03
PDS-1RD110-20	p.03
PDS-1RD110-30	p.03
PDS-1RD110-40	p.03
PDS-1RD130-20	p.03
PDS-1RD130-30	p.03
PDS-1RD130-40	p.03
PDS-1RD170-35	p.03
PDS-1RD170-40	p.03
PDS-1RD170-50	p.03
PDS-1RD240-30	p.03
PDS-1RD240-40	p.03
PDS-1RD240-45	p.03
PDS-1RD240-50	p.03

Order Code	Page
PDS-1RF110-30	p.04
PDS-1RF110-40	p.04
PDS-1RF130-30	p.04
PDS-1RF130-40	p.04
PDS-1RF170-35	p.04
PDS-1RF170-40	p.04
PDS-1RF240-40	p.04
PDS-1RTCD110-45	p.03
PDS-1RTCD130-45	p.03
PDS-1RTD110-45	p.03
PDS-1RTD130-45	p.03
PDS-2AES-50	p.10
PDS-2AES-60	p.10
PDS-2AES-75	p.10
PDS-2AMS-50	p.10
PDS-2AMS-60	p.10
PDS-2AMS-75	p.10
PDS-2AS-50	p.10
PDS-2AS-60	p.10
PDS-2AS-75	p.10
PDS-2AUL-50	p.06
PDS-2AUL-60	p.06
PDS-2AUL-75	p.06
PDS-2BCDL-30	p.09
PDS-2BCDL-30E	p.09
PDS-2BCDL-50	p.09
PDS-2BCDL-50E	p.09
PDS-2BCDUL-30	p.05
PDS-2BCDUL-50	p.05
PDS-2BES-40	p.10
PDS-2BES-60	p.10
PDS-2BMS-40	p.10
PDS-2BMS-60	p.10
PDS-2BS-40	p.10
PDS-2BS-60	p.10
PDS-2CBES-20	p.09
PDS-2CBES-30	p.09
PDS-2CBES-40	p.09
PDS-2CBES-50	p.09
PDS-2CBES-60	p.09
PDS-2CBES-70	p.09
PDS-2CBL-20	p.09
PDS-2CBL-30	p.09
PDS-2CBL-40	p.09
PDS-2CBL-50	p.09
PDS-2CBL-60	p.09
PDS-2CBL-70	p.09
PDS-2CBM-20	p.09
PDS-2CBM-30	p.09
PDS-2CBM-40	p.09

Order Code	Page
PDS-2CBM-50	p.09
PDS-2CBM-60	p.09
PDS-2CBM-70	p.09
PDS-2CBMS-20	p.09
PDS-2CBMS-30	p.09
PDS-2CBMS-40	p.09
PDS-2CBMS-50	p.09
PDS-2CBMS-60	p.09
PDS-2CBMS-70	p.09
PDS-2CBS-20	p.09
PDS-2CBS-30	p.09
PDS-2CBS-40	p.09
PDS-2CBS-50	p.09
PDS-2CBS-60	p.09
PDS-2CBS-70	p.09
PDS-2CDES-10	p.08
PDS-2CDES-10E	p.08
PDS-2CDES-15	p.08
PDS-2CDES-15E	p.08
PDS-2CDES-20	p.08
PDS-2CDES-30	p.08
PDS-2CDES-40	p.08
PDS-2CDES-50	p.08
PDS-2CDES-60	p.08
PDS-2CDES-70	p.08
PDS-2CDL-20E	p.08
PDS-2CDL-30	p.08
PDS-2CDL-30E	p.08
PDS-2CDL-40	p.08
PDS-2CDL-40E	p.08
PDS-2CDL-50	p.08
PDS-2CDL-50E	p.08
PDS-2CDL-70	p.08
PDS-2CDM-20	p.08
PDS-2CDM-20E	p.08
PDS-2CDM-30	p.08
PDS-2CDM-30E	p.08
PDS-2CDM-40	p.08
PDS-2CDM-40E	p.08
PDS-2CDM-50	p.08
PDS-2CDM-50E	p.08
PDS-2CDM-60	p.08
PDS-2CDM-70	p.08
PDS-2CDMS-10	p.08
PDS-2CDMS-15	p.08
PDS-2CDMS-20	p.08
PDS-2CDMS-30	p.08
PDS-2CDMS-40	p.08

Order Code	Page
PDS-2CDMS-50	p.08
PDS-2CDMS-60	p.08
PDS-2CDMS-70	p.08
PDS-2CDS-10	p.08
PDS-2CDS-15	p.08
PDS-2CDS-20	p.08
PDS-2CDS-20E	p.08
PDS-2CDS-30	p.08
PDS-2CDS-30E	p.08
PDS-2CDS-40	p.08
PDS-2CDS-40E	p.08
PDS-2CDS-50	p.08
PDS-2CDS-50E	p.08
PDS-2CDS-60	p.08
PDS-2CDS-70	p.08
PDS-2CES-10	p.09
PDS-2CES-20	p.09
PDS-2CES-30	p.09
PDS-2CES-40	p.09
PDS-2CES-50	p.09
PDS-2CES-60	p.09
PDS-2CES-70	p.09
PDS-2CL-20	p.09
PDS-2CL-30	p.09
PDS-2CL-40	p.09
PDS-2CL-50	p.09
PDS-2CL-60	p.09
PDS-2CL-70	p.09
PDS-2CM-10	p.09
PDS-2CM-20	p.09
PDS-2CM-30	p.09
PDS-2CM-40	p.09
PDS-2CM-50	p.09
PDS-2CM-60	p.09
PDS-2CM-70	p.09
PDS-2CMS-10	p.09
PDS-2CMS-20	p.09
PDS-2CMS-30	p.09
PDS-2CMS-40	p.09
PDS-2CMS-50	p.09
PDS-2CMS-60	p.09
PDS-2CMS-70	p.09
PDS-2CS-10	p.09
PDS-2CS-20	p.09
PDS-2CS-30	p.09
PDS-2CS-40	p.09
PDS-2CS-50	p.09
PDS-2CS-60	p.09
PDS-2CS-70	p.09
PDS-2CYL-50	p.10

Order Code	Page
PDS-2CYL-60	p.10
PDS-2CYUL-50	p.06
PDS-2CYUL-60	p.06
PDS-2DES-06	p.08
PDS-2DES-10	p.08
PDS-2DES-15	p.08
PDS-2DES-20	p.08
PDS-2DES-30	p.08
PDS-2DES-40	p.08
PDS-2DES-50	p.08
PDS-2DES-60	p.08
PDS-2DES-70	p.08
PDS-2DL-20	p.08
PDS-2DL-30	p.08
PDS-2DL-40	p.08
PDS-2DL-50	p.08
PDS-2DL-60	p.08
PDS-2DL-70	p.08
PDS-2DM-10	p.08
PDS-2DM-20	p.08
PDS-2DM-30	p.08
PDS-2DM-40	p.08
PDS-2DM-50	p.08
PDS-2DM-60	p.08
PDS-2DM-70	p.08
PDS-2DMS-06	p.08
PDS-2DMS-10	p.08
PDS-2DMS-15	p.08
PDS-2DMS-20	p.08
PDS-2DMS-30	p.08
PDS-2DMS-40	p.08
PDS-2DMS-50	p.08
PDS-2DMS-60	p.08
PDS-2DMS-70	p.08
PDS-2DS-06	p.08
PDS-2DS-10	p.08
PDS-2DS-15	p.08
PDS-2DS-20	p.08
PDS-2DS-30	p.08
PDS-2DS-40	p.08
PDS-2DS-50	p.08
PDS-2DS-60	p.08
PDS-2DS-70	p.08
PDS-2ECDES-30	p.08
PDS-2ECDES-40	p.08
PDS-2ECDES-50	p.08
PDS-2ECDES-60	p.08
PDS-2ECDL-30	p.08
PDS-2ECDL-40	p.08
PDS-2ECDL-50	p.08

Order Code	Page
PDS-2ECDL-60	p.08
PDS-2ECDM-30	p.08
PDS-2ECDM-40	p.08
PDS-2ECDM-50	p.08
PDS-2ECDM-60	p.08
PDS-2ECDMS-30	p.08
PDS-2ECDMS-40	p.08
PDS-2ECDMS-50	p.08
PDS-2ECDMS-60	p.08
PDS-2ECDMS-70	p.08
PDS-2ECDS-30	p.08
PDS-2ECDS-40	p.08
PDS-2ECDS-50	p.08
PDS-2ECDS-60	p.08
PDS-2EES-40	p.10
PDS-2EES-60	p.10
PDS-2EFUL-40	p.06
PDS-2EFUL-60	p.06
PDS-2EL-40	p.10
PDS-2EL-60	p.10
PDS-2EMS-40	p.10
PDS-2EMS-60	p.10
PDS-2ES-40	p.10
PDS-2ES-60	p.10
PDS-2FES-12	p.10
PDS-2FES-16	p.10
PDS-2FES-21	p.10
PDS-2FMS-12	p.10
PDS-2FMS-16	p.10
PDS-2FMS-21	p.10
PDS-2FS-12	p.10
PDS-2FS-16	p.10
PDS-2FS-21	p.10
PDS-2LES-11	p.10
PDS-2LES-26	p.10
PDS-2LMS-11	p.10
PDS-2LMS-26	p.10
PDS-2LS-11	p.10
PDS-2LS-26	p.10
PDS-2MHCDUL-20	p.05
PDS-2MHCDUL-30	p.05
PDS-2MHCDUL-40	p.05
PDS-2MHFUL-18	p.06
PDS-2MHFUL-22	p.06
PDS-2MHFUL-30	p.06
PDS-2NCBES-18	p.10
PDS-2NCBES-22	p.10
PDS-2NCBES-30	p.10
PDS-2NCBL-18	p.10
PDS-2NCBL-22	p.10

Order Code	Page
PDS-2NCBL-30	p.10
PDS-2NCBM-18	p.10
PDS-2NCBM-22	p.10
PDS-2NCBM-30	p.10
PDS-2NCBMS-18	p.10
PDS-2NCBMS-22	p.10
PDS-2NCBMS-30	p.10
PDS-2NCBS-18	p.10
PDS-2NCBS-22	p.10
PDS-2NCBS-30	p.10
PDS-2NES-18	p.09
PDS-2NES-22	p.09
PDS-2NES-30	p.09
PDS-2NL-18	p.09
PDS-2NL-22	p.09
PDS-2NL-30	p.09
PDS-2NM-18	p.09
PDS-2NM-22	p.09
PDS-2NM-30	p.09
PDS-2NMS-18	p.09
PDS-2NMS-22	p.09
PDS-2NMS-30	p.09
PDS-2NS-18	p.09
PDS-2NS-22	p.09
PDS-2NS-30	p.09
PDS-2PES-60	p.10
PDS-2PMS-60	p.10
PDS-2PS-60	p.10
PDS-2QEUL-30	p.06
PDS-2QEUL-40	p.06
PDS-2QEUL-50	p.06
PDS-2QEUL-60	p.06
PDS-2RCDUL-10	p.05
PDS-2RCDUL-15	p.05
PDS-2RCDUL-20	p.05
PDS-2RCDUL-25	p.05
PDS-2RCDUL-30	p.05
PDS-2RCDUL-40	p.05
PDS-2RCDUL-50	p.05
PDS-2RCDUL-60	p.05
PDS-2RCDUL-70	p.05
PDS-2RDUL-06	p.05
PDS-2RDUL-20	p.05
PDS-2RDUL-25	p.05
PDS-2RDUL-30	p.05
PDS-2RDUL-35	p.05
PDS-2RDUL-40	p.05
PDS-2RDUL-50	p.05
PDS-2RDUL-60	p.05
PDS-2RECDUL-30	p.05

Index

Order Code	Page
PDS-2RECDUL-40	p.05
PDS-2RECDUL-50	p.05
PDS-2RECDUL-60	p.05
PDS-2RECDUL-70	p.05
PDS-2RES-13	p.10
PDS-2RFUL-10	p.06
PDS-2RFUL-15	p.06
PDS-2RFUL-20	p.06
PDS-2RFUL-25	p.06
PDS-2RFUL-30	p.06
PDS-2RFUL-35	p.06
PDS-2RFUL-40	p.06
PDS-2RFUL-50	p.06
PDS-2RFUL-60	p.06
PDS-2RFUL-70	p.06
PDS-2RFUL-88	p.06
PDS-2RMS-13	p.10
PDS-2RSCDUL-30	p.05
PDS-2RSCDUL-40	p.05
PDS-2RSCDUL-50	p.05
PDS-2RSCDUL-60	p.05
PDS-2RUL-13	p.06
PDS-2SCDES-30	p.09
PDS-2SCDES-40	p.09
PDS-2SCDES-50	p.09
PDS-2SCDES-60	p.09
PDS-2SCDL-30	p.09
PDS-2SCDL-40	p.09
PDS-2SCDL-50	p.09
PDS-2SCDL-60	p.09
PDS-2SCDM-30	p.09
PDS-2SCDM-40	p.09
PDS-2SCDM-50	p.09
PDS-2SCDM-60	p.09
PDS-2SCDMS-30	p.09
PDS-2SCDMS-40	p.09
PDS-2SCDMS-50	p.09
PDS-2SCDMS-60	p.09
PDS-2SCDS-30	p.09
PDS-2SCDS-40	p.09
PDS-2SCDS-50	p.09
PDS-2SCDS-60	p.09
PDS-2TD-10	p.11
PDS-2TD-15	p.11
PDS-2TD-20	p.11
PDS-2WES-10	p.10
PDS-2WES-12	p.10
PDS-2WES-15	p.10
PDS-2WMS-10	p.10
PDS-2WMS-12	p.10

Order Code	Page
PDS-2WMS-15	p.10
PDS-2WS-10	p.10
PDS-2WS-12	p.10
PDS-2WS-15	p.10
PDS-2WUL-10	p.06
PDS-2WUL-12	p.06
PDS-2WUL-15	p.06
PDS-3CDL-20	p.11
PDS-3CDL-30	p.11
PDS-3CDL-30E	p.11
PDS-3CDL-40	p.11
PDS-3CDL-40E	p.11
PDS-3CDL-50	p.11
PDS-3CDL-50E	p.11
PDS-3CDL-60	p.11
PDS-3CDL-70	p.11
PDS-3CDM-20	p.11
PDS-3CDM-30	p.11
PDS-3CDM-30E	p.11
PDS-3CDM-40	p.11
PDS-3CDM-40E	p.11
PDS-3CDM-50	p.11
PDS-3CDM-50E	p.11
PDS-3CDM-60	p.11
PDS-3CDM-70	p.11
PDS-3CDS-20	p.11
PDS-3CDS-30	p.11
PDS-3CDS-40	p.11
PDS-3CDS-50	p.11
PDS-3CDS-60	p.11
PDS-3CDS-70	p.11
PDS-3CL-20	p.11
PDS-3CL-30	p.11
PDS-3CL-40	p.11
PDS-3CL-50	p.11
PDS-3CL-60	p.11
PDS-3CL-70	p.11
PDS-3CM-10	p.11
PDS-3CM-20	p.11
PDS-3CM-30	p.11
PDS-3CM-40	p.11
PDS-3CM-50	p.11
PDS-3CM-60	p.11
PDS-3CM-70	p.11
PDS-3CS-10	p.11
PDS-3CS-20	p.11
PDS-3CS-30	p.11
PDS-3CS-40	p.11
PDS-3CS-50	p.11
PDS-3CS-60	p.11

Order Code	Page
PDS-3CS-70	p.11
PDS-3DL-20	p.11
PDS-3DL-20E	p.11
PDS-3DL-30	p.11
PDS-3DL-40	p.11
PDS-3DL-50	p.11
PDS-3DL-60	p.11
PDS-3DL-70	p.11
PDS-3DM-10	p.11
PDS-3DM-20	p.11
PDS-3DM-20E	p.11
PDS-3DM-30	p.11
PDS-3DM-40	p.11
PDS-3DM-50	p.11
PDS-3DM-60	p.11
PDS-3DM-70	p.11
PDS-3DS-10	p.11
PDS-3DS-20	p.11
PDS-3DS-30	p.11
PDS-3DS-40	p.11
PDS-3DS-50	p.11
PDS-3DS-60	p.11
PDS-3DS-70	p.11
PDS-3MCS-CB16	p.12
PDS-3MCS-CB30	p.12
PDS-3MCS-DD254	p.12
PDS-3NL-18	p.11
PDS-3NL-22	p.11
PDS-3NL-30	p.11
PDS-3NM-18	p.11
PDS-3NM-22	p.11
PDS-3NM-30	p.11
PDS-3NS-18	p.11
PDS-3NS-22	p.11
PDS-3NS-30	p.11
PDS-CRA-L	p.12
PDS-CRA-M	p.12
PDS-CRA-P	p.12
PDS-OSC-422	p.16
PDS-OSC-432	p.16
PDS-RE-4023	p.13
PDS-RE-4023E	p.13
PDS-RE-4028	p.13
PDS-RE-4028E	p.13
PDS-RE-4114	p.13
PDS-RE-4127	p.13
PDS-RE-4217	p.13
PDS-RE-4223	p.13
PDS-RE-6127	p.13
PDS-RE-R011	p.14

Order Code	Page
PDS-RE-R012	p.14
PDS-RE-R111	p.14
PDS-RE-R112	p.14
PDS-ROU-L	p.12
PDS-ROU-M	p.12
PDS-ROU-P	p.12
PDS-SO-411	p.15
PDS-SO-412	p.15
PDS-SO-422	p.15
PDS-SO-431	p.15
PDS-SO-432	p.15
PDS-SO-433	p.15
PDS-SO-441	p.15
PDS-SO-442	p.15
PDS-SO-452	p.15
PDS-SO-453	p.15
PDS-SO-632	p.15
PDS-SO-642	p.15
PDS-SO-653	p.15
PDS-WSP-L	p.12
PDS-WSP-M	p.12

